

LABORATORY SERVICES BUREAU		
Document: Toxicology Procedures	Policy Number: 1265	Revision: 13
Subject: TOX-SOP-37 Protocol for the Analysis of Phencyclidine in Blood	Approved: Gallegos, Amanda	
PHOENIX POLICE DEPARTMENT	Effective: 3/11/2024 5:19:26 PM	Page 1 of 4

1. PROTOCOL FOR THE ANALYSIS OF PHENCYCLIDINE IN BLOOD

PURPOSE

The following method describes the quantitation of phencyclidine in blood, serum, plasma or other biological samples by GC/MS. Samples which have screened positive by a preliminary test, as well as special requests or retest requests will follow the following protocol.

PLAN

A. Equipment:

- (1) GC/MS with a 5% diphenylpolysiloxane, 95% dimethylpolysiloxane, 15/30 meter, 0.25 micron film thickness column
- (2) Positive Pressure Manifold
- (3) SPE Column – Silica gel - Co polymeric bonded phase with a hydrophobic cation exchange (CSDAU203 or XRDAH203)
- (4) Sample concentrator with UHP nitrogen
- (5) Centrifuge
- (6) Vortex Mixer / Multi-tube vortex mixer

B. Reagents:

- (1) **100 mM Phosphate buffer solution.** Dissolve 1.70 grams of Na₂HPO₄ and 12.14 g NaH₂PO₄·H₂O in 800 ml of deionized water. Dilute to 1000 ml with deionized water. Mix well. pH should be 5.5-6.0. If necessary, adjust with 100 mM monobasic sodium phosphate (lowers pH) or 100 mM dibasic sodium phosphate (raises pH). Store refrigerated. Stable for six months.
- (2) **Methanol.** Prepare a transfer bottle of ACS/HPLC grade methanol. Label accordingly. Store in glass at room temperature. Stable until consumed.
- (3) **100 mM Hydrochloric Acid (HCl).** To 400 ml of deionized water, add 8.4 ml concentrated HCl. Dilute to 1 L with deionized water. Mix well. Store at room temperature. Stable for 2 years.
- (4) **78:20:2 methylene chloride: isopropanol: ammonium hydroxide Elution Solvent.** Prepare fresh daily. Add ammonium hydroxide to isopropanol, followed by methylene chloride (i.e. per 10 mL of elution solvent add approximately 200 µL ammonium hydroxide). Mix thoroughly, elution solvent should have a turbid appearance when thoroughly mixed.
- (5) **Ethyl acetate.** Prepare a transfer bottle of ACS/HPLC grade ethyl acetate. Label accordingly. Store in glass at room temperature. Stable until consumed.
- (6) **Deionized Water** (DI Water) Label, stable until consumed.

LABORATORY SERVICES BUREAU

Document: Toxicology Procedures

Policy Number:
1265

Revision:
13

Subject: TOX-SOP-37 Protocol for the Analysis of Phencyclidine in Blood

Approved:
Gallegos, Amanda

PHOENIX POLICE DEPARTMENT

Effective: 3/11/2024 5:19:26 PM

Page 2 of 4

C. Standards (Store refrigerated. Stable 2 years if prepared in house or per manufacturer's recommendation):

- (1) **1 mg/ml Phencyclidine HCl stock standard.** Prepare by weighing out 11.5 mg phencyclidine HCl and dissolving in 10 ml of methanol to arrive at a concentration of 1.0 mg/ml. Or purchase 1.0 mg/ml ampoule (Cerilliant P-007).
- (2) **100 µg/ml stock D5-phencyclidine internal standard.** Purchase a 100 µg/ml ampoule (Cerilliant P-003).

D. Calibrators and Internal Standard (Store refrigerated. Stable for 2 years):

- (1) **1 ng/µl Phencyclidine calibrator stock solution in methanol.** To a 10 ml volumetric flask add 10 µl of 1 mg/ml stock PCP standard. Dilute to volume with methanol
- (2) **1 ng/µl D5-PCP internal standard.** In a 10 ml volumetric flask add 100 µl of 100 µg/ml D5-PCP stock standard. Dilute to volume with methanol. (May also prepare larger stock volume by adjusting volumes accordingly to account for equivalent concentrations of each analyte.)

E. Quality Controls: (Store refrigerated)

- (1) **Positive Controls.** Prepared on day of use from a different lot of stock solution than that used to prepare calibrators or purchased from an external vendor. Target concentration is 30, 60, and 80 ng/ml for PCP. Additional controls shall be prepared when appropriate, to coincide with any limited sample volumes and/or dilution of case samples.
- (2) **Negative Control.** Blank blood prepared in house consisting of 50% saline, 50% packed red blood cells, and 5g sodium fluoride/1g potassium oxalate (per 500 ml prepared blood) will be used as negative control.

F. Solid Phase Extraction (SPE):

(1) **Sample Preparation.**

Prepare in appropriately labeled culture tubes as follows:

Prepare a set of calibrators at 10, 25, 50, and 100ng/ml using the above calibrator stocks in 1 mL of blank blood along with 1mL* negative control, positive controls, case samples and additional controls at appropriate volumes and/or dilutions if applicable. Add 50µl* of working internal standard to each tube.

*In case samples where a limited amount of blood is received use the same fraction of internal standard as the sample, as an example ½ ml of blood and 25 µl of internal standard. Case samples which must be diluted to fall within the calibration range will receive full internal standard, as an example x2 by using ½ mL sample with 50 µl of internal standard.

- (a) Add 1.5 ml of water and vortex each tube until thoroughly mixed.
- (b) Add 1.0 ml of 100 mM phosphate buffer and vortex until thoroughly mixed.
- (c) Centrifuge at 3,500 rpm for 5 minutes.

LABORATORY SERVICES BUREAU

Document: Toxicology Procedures

Policy Number:
1265

Revision:
13

Subject: TOX-SOP-37 Protocol for the Analysis of Phencyclidine in Blood

Approved:
Gallegos, Amanda

PHOENIX POLICE DEPARTMENT

Effective: 3/11/2024 5:19:26 PM

Page 3 of 4

(2) **Column Conditioning**

Pass through the column sequentially the following reagents at <1.0 ml/min or by gravity only:

- (a) 2 ml of methanol
- (b) 2 ml of deionized water
- (c) 1 ml of 100 mM phosphate buffer

Take care to prevent sorbent from drying out.

(3) **Sample Application**

Apply sample to column, being careful to not allow the sediment, if present, which will be in the base of the centrifuge tube to pass. Flow rate should be about 1.0 ml/minute or gravity only.

(4) **Column Rinse and Elution**

Pass through the column sequentially the following reagents, at 1-2 ml/min:

- (a) 3 ml of deionized water.
- (b) 1 ml of 100 mM HCl
- (c) 2 ml methanol
- (d) Dry column under full pressure (20-25psi) for 15 minutes.
- (e) Elute by gravity with 1.5 ml of freshly prepared 78:20:2 methylene chloride: isopropanol: ammonium hydroxide solution directly into appropriately labeled microvials.
- (f) Evaporate to dryness. Reconstitute with 60 µl of ethyl acetate, vortex (transfer to vial insert if using screw cap autosampler vials) and cap.

G. Data Acquisition and Analysis:

- (1) Make sure the Autotune was performed, rinse vials filled, etc
- (2) Set up a sequence with the calibrators injected first in order to calibrate the instrument used. Subsequent injections to include negative control, positive controls (positive controls should be included throughout the batch, i.e. beginning, mid-run and end of run when possible) and solvent blanks prior to case samples. For samples requiring dilution, add the appropriate sample multiplier in the sequence table. Load samples onto autosampler according to sequence and have it verified by another analyst before or after analysis but prior to unloading.
- (3) The ion ratios and retention times should be set by a mid-level calibrator.
- (4) Analyze using the appropriate method on GC/MS.

H. Results and Acceptability:

- (1) Calibration $R^2 \geq 0.99$ and calibrators within 20% of set value.

LABORATORY SERVICES BUREAU		
Document: Toxicology Procedures	Policy Number: 1265	Revision: 13
Subject: TOX-SOP-37 Protocol for the Analysis of Phencyclidine in Blood	Approved: Gallegos, Amanda	
PHOENIX POLICE DEPARTMENT	Effective: 3/11/2024 5:19:26 PM	Page 4 of 4

- (2) Positive control within 20% of target concentration.
- (3) If the above two criteria are not met the analyte may be reported qualitatively
- (4) Negative control < 25% of area count of cutoff calibrator
- (5) Qualifier ion ratios within 20% as set from calibrator.
- (6) Retention time within 2% as set from calibrator.
- (7) Chromatographically acceptable i.e. peak purity ≥90% for target/quantitative ion.
- (8) Blank prior to sample < 25% area count of cutoff calibrator
- (9) Quantitation ≥ 10 ng/ml; results greater than highest calibrator will be reported qualitatively as such, and samples which included a dilution factor will be reported greater than the highest calibrator multiplied by the applicable dilution factor.
- (10) Results will be truncated and documented in case notes to two significant figures.