

LABORATORY SERVICES BUREAU

Document: Controlled Substances Analysis Manual

Policy Number:
1520

Revision:
6

Subject: CS-SOP-40 Morphine


Approved:
Bell, Erica

PHOENIX POLICE DEPARTMENT Effective: 7/15/2024 2:23:29 PM

Page 1 of 3

1. MORPHINE

A. Structure, Empirical Formula, Molecular Weight



Morphine
 $C_{17}H_{19}NO_3$
MW 285.3

B. Synonyms: NA

C. Trade Names: Duramorph®, MS Contin®, MSIR®, Oramorph SR™, Roxanol™, Kadian®, etc.

D. Drug Action: Narcotic analgesic

E. Common pharmaceutical/street forms: 15-200 mg tablets, 15-100 mg capsules, 2 to 20 mg/mL oral solution, 0.5 to 25 mg/mL injectable solution, 5 to 30 mg rectal suppositories. Morphine is present in opium in concentrations of 9-15%.

F. Solubility:

- (1) Free base: Slightly soluble in isopropyl alcohol, ethanol, methanol, and chloroform, insoluble in water, and diethyl ether
- (2) Various salts: Water, slightly soluble in ethanol and methanol, insoluble in chloroform and diethyl ether

G. Extraction:

- (1) Dry extraction
 - (a) Determine amount of sample required for extraction based on the concentration of morphine in the sample.
 - (b) Grind sample and place in spot plate or test tube.
 - (c) Extract with methanol:isopropyl alcohol (4:1).
 - (d) Filter to remove any insoluble material if necessary.
 - (e) Evaporate the extraction solvent and reconstitute in C15 methanol.

LABORATORY SERVICES BUREAU		
Document: Controlled Substances Analysis Manual	Policy Number: 1520	Revision: 6
Subject: CS-SOP-40 Morphine	Approved: Bell, Erica	
PHOENIX POLICE DEPARTMENT	Effective: 7/15/2024 2:23:29 PM	Page 2 of 3

(2) Dry extraction

- (a) Determine amount of sample required for extraction based on the concentration of morphine in the sample.
- (b) Grind sample and place in spot plate or test tube.
- (c) Extract with C15 methanol.
- (d) Filter to remove any insoluble material if necessary.

(3) Dry extraction for controlled release tablets

- (a) Determine amount of sample required for extraction based on the concentration of morphine in the sample.
- (b) Ether wash the sample and discard the ether.
- (c) Extract the remaining sample with C15 methanol.
- (d) Filter to remove any insoluble material if necessary.

(4) Acid/base extraction for samples combined with acetaminophen or aspirin

- (a) Determine amount of sample required for extraction based on the concentration of morphine in the sample.
- (b) Grind tablet or remove liquid and place in a test tube.
- (c) Add 0.2 N sulfuric acid.
- (d) Extract with chloroform two times, discarding the chloroform (bottom layer). If emulsions occur, the sample may be centrifuged.
- (e) Basify (saturate) the acidic solution with solid potassium carbonate. Verify that the pH is basic.
- (f) Extract with chloroform:isopropyl alcohol (4:1). If emulsions occur, the sample may be centrifuged.
- (g) Remove the chloroform:isopropyl alcohol layer and place into a clean well.
- (h) Repeat steps f and g, combining the chloroform:isopropyl alcohol extracts.

LABORATORY SERVICES BUREAU		
Document: Controlled Substances Analysis Manual	Policy Number: 1520	Revision: 6
Subject: CS-SOP-40 Morphine	Approved: Bell, Erica	
PHOENIX POLICE DEPARTMENT	Effective: 7/15/2024 2:23:29 PM	Page 3 of 3

H. Chemical indicator tests:

- (1) Marquis: Purple
- (2) Liebermann's: Brown-black
- (3) Froehde's: Purple
- (4) Mecke's: Green

I. TLC:

- (1) Mobile Phase:
 - (a) System 1: Chloroform:methanol (90:10)
 - (b) System 2: Methanol:conc. ammonium hydroxide (100:1.5)
- (2) Locator: Acidified iodoplatinate, Dragendorff reagent, Marquis reagent

J. GC/MS: Analyze using "Drugs1" program.

K. Comments: Used for pain relief.

L. Report as: Morphine, a narcotic drug.

M. References:

- (1) Analytical Profiles of Narcotic Analgesics, CND Analytical, Inc., Auburn, AL, 1991.